

Connected to

... REPRODUCTION BIOTECHNOLOGIES

α pha

BETTER INSEMINATION

- High Security Insemination
- Animal Well-being
- Time Saving

Alpha



BETTER INSEMINATION

High Security Insemination

- Palpable dome head design, facilitates optimum placement
- Better straw fit, no backflow
- Crystalline plastic, easy reading of straw print including through the head
- Dual-lateral delivery, improves hygiene
- Transparent head, makes fluids easy to see upon withdrawal
- Improved connection, prevents separation from the gun
- Flexible and resilient, assists in placement

Animal Well-being

- Polished dome head, passes smoothly into cervix and through cervical folds
- Ultra-smooth, atraumatic

Time Saving

- Insert integrated in the head, quick loading
- Two lateral exit holes, for reliable and consistent deposit of semen
- Crystalline plastic, allows barcode scanning
- Unsplit, connects to Kombicolor or spiral gun
- Universal, fits mini or medium straws

Optimized for Kombicolor

- Inseminate 4 cm deeper than other guns
- Fast connection and easy removal of the sheath
- Quick identification with 5 different colored rings
- Universal, mini and medium straws
- Non-magnetic, austenitic 304L 18/10 stainless steel for years of use
- Easy to clean for improved hygiene



Ref. 024485

Packaged in bag of 50. Master box of 1500

