



www.imv-technologies.com

# EasyClean™

## A Washing Solution for Flow Cytometer Specific for Semen Analysis

### A Combination of adapted agents

- Increases service life of capillary
- Can be used within regular washing procedure
- Preventive and curative treatment
- Ready-to-use

### Specific for semen analysis

- Removes any residue of proteins, cells or fluorochromes
- Formula developed for egg yolk-based media and clear media
- Prevents and eliminates cloggings
- Proven effectiveness

### Optimized bottle design

- Dropper to fill tubes and washing bottle easily
- Eliminates pipetting process: reduces risks of solution contamination

100 ml Bottle  
Packaged in lots of 4  
Ref. 024760



- ISO 9001: 2008 Certification and ISO 13 485: 2003 Medical Standard
- A production site complying with the most rigorous quality and safety standards
- Continuous improvement of products, quality and control processes
- All materials and finished products potentially in contact with living cells are tested for biocompatibility
- 95% made in France



**imv**  
TECHNOLOGIES



+ 33 (0) 233 346 464  
contact@imv-technologies.com

www.imv-technologies.com



communication 11/2012 - Non contractual pictures and specifications