

Connected to

... REPRODUCTION BIOTECHNOLOGIES

GTB 250

MANUAL BAG FILLING AND SEALING MACHINE

Our entry level packaging system

- **Efficient**
Fills 6 bags at a time
250 bags per hour
- **Simple and user-friendly**
Fills by gravity
Easy to seal
- **Perfect sanitation**
Uses disposable tubing and filling needles



GTB 250 specifications

References
Compatible bags
Throughput (Bag/H)
Sealing
Filling
Bag advance
Bag cutting
Labelling
Footprint L x W x H
Material
Electric connection
Power
Weight

GTB 250

022267 (230V) / 022322 (110V)
GTB bag (021595) / Bactibag 250 (027146)
250
Impulse Sealer
Manual
Manual
Manual
Manual
100x40x65cm / 39"x16"x25"
Stainless steel and aluminum
230/110 – 50/60 Hz
150 W
16.8 kg / 37 lbs

GTB 250 spare parts and consumables

References	Description	Quantity
022393	Sealing Wire	2
015404	Teflon Strip	2
007730	Sinker Weight	1
006447	Filling nozzle & tubing (bag of 5)	1
005614	Filling clamp	2
023265	Priming bottles	5
021595	GTB bags	500

