



Kombicolor: when size matters

4 cm longer than conventional AI guns*



Only 2 parts

Sheath easily secured

Bigger ring fits better your fingers


One gun for mini & medium straws

5 colors available for better identification

Stainless steel built to last


*with sheath

WORLD LEADER IN REPRODUCTION BIOTECHNOLOGIES



www.imv-technologies.com

Kombicolor



STATE OF THE ART INNOVATION

No more loose or lost rings
 A specific part secures the sheath allowing a very easy assembling and disassembling
 Best 18/10 stainless steel available for years of trouble-free use
 Easy to clean for a perfect hygiene

HOW DOES IT WORK?

Thanks to a new fixation, the entire length of the sheath is optimized for deeper insemination

Fingers are securely kept away on a bigger colored plastic flange

To be used with IMV unsplit sheath



WHAT ARE THE BENEFITS?

Time saving: very quick sheath securing & removal.

One gun fits mini & medium straws

5 colors available for clear & quick identification

